

Resuscitation



OxyAL™ Lifesaver

User-friendly resuscitation kit for swift response in an emergency

OxyAL™ Lifesaver is a compact resuscitation kit especially designed by Air Liquide for quick emergency response.

- ✓ Compact and easy to carry
- ✓ Handy storage for accessories
- ✓ Available with reusable or disposable resuscitator kit

Item Description	Item Number
OxyAL with Reusable resuscitator kit (Reusable adult resuscitator with 1500ml ventilation bag, 60cmH2O pop-off valve, size 5 silicone reusable mask, 2500ml reservoir bag & 2.1m non-kink tubing)	3463SK
OxyAL with Disposable resuscitator kit (Disposable adult resuscitator with 1500ml ventilation bag, 60cmH2O pop-off valve, size 5 disposable mask, 2500ml reservoir bag & 2.1m non-kink tubing)	3468SK



Parts & Accessories	Item Number
Disposable Resuscitator	2858S
Reusable Resuscitator	2857S
OxyAL Carry Bag	1903S
Guedel airways; Size 0 (60mm)	1893S
3 (90mm)	1894S
5 (110mm)	1895S
Oxygen tube, 7 ft	843S
O2 MASK Adult	154S
Paediatric	157S

Air Liquide Australia

🌐 au.airliquide.com

✉ auqueries@airliquide.com

📞 For queries: 1300 733 011

📞 For emergencies: 1800 812 588